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
# Learn.Genetics (Genetic Science Learning Center)

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<http://learn.genetics.utah.edu/content/labs/extraction/>

### VIRTUAL LAB: DNA EXTRACTION

DNA is extracted from human cells for a variety of reasons. With a pure sample of DNA you can test a newborn for a genetic disease, analyze forensic evidence, or study a gene involved in cancer.

Try this virtual laboratory to perform a cheek swab and extract DNA from human cells.

➤ Why, you may ask, do you need, human DNA?

.....  
.....  
.....

➤ The steps you will follow to purify DNA from a cheek swab are:

- 1).....
- 2).....
- 3).....
- 4).....

➤ Put in the following table the materials and equipment you will need to purify DNA from cheek cells:

materials	equipment

➤ Build up a glossary by matching words and their translation:

words	translation	words	translation

*Basic Handout*

The steps to extract the DNA are as follows:

- 1) Grind the cheek cells
- 2) Heat cells open to release DNA
- 3) Separate DNA from proteins and other
- 4) Isolate denatured DNA

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1) Grind the cheek cells

2) Heat cells open to release DNA

3) Separate DNA from proteins and other

4) Isolate denatured DNA

5) Grind the cheek cells

6) Heat cells open to release DNA

7) Separate DNA from proteins and other

8) Isolate denatured DNA

9) Grind the cheek cells

10) Heat cells open to release DNA

11) Separate DNA from proteins and other

12) Isolate denatured DNA

Conclusion

This experience has been very useful in the practical and technological field to understand the extraction of DNA in an effective way.

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Conclusion

This experience has been very useful because with the DNA we can see how it is made and how it is used. It is also very interesting to see how the DNA is made and how it is used. In English we have learned a lot of knowledge of the English language.

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**LYSIS SOLUTION**  
 (Separate of DNA from histones)

**CONCENTRATED SALT SOLUTIONS**  
 (soluzione salina)  
 precipitates proteins and debris.

**ISOPROPYL ALCOHOL**  
 (alcohol isopropilico)  
 Makes it not soluble DNA

**BUCCAL SWAB**  
 (tampona buccale)  
 it takes cells from skin

**MICROPIPETTES**  
 (micropipette)  
 Small pipette that picks up DNA

It encountered any problems.  
 has helped us to learn English

Handwritten notes: "is written in our class... in second moment... we have extracted the DNA... it is quite convertible DNA"

Diagrams: A pipette, a microcentrifuge tube, and a microcentrifuge.

**BUCCAL SWAB**  
 (to extract cells)

**MICROCENTRIFUGE**  
 (contains material)

**WATER BATH**  
 (HEATS THE SOLUTIONS)

**MICROPIPETTES**  
 (TAKE SOLUTION)

**TRIS**

**CENTRIFUGE**  
 (SEPARATE SOLUTIONS)

**ETHANOL**  
 (SOLUTION)

**LYSIS SOLUTION**  
 (SEPARATE DNA FROM THE CELLULAR MEMBRANE)

**ISOPROPYL ALCOHOL**  
 (PRECIPITATION OF DNA)

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**Crime Scene Investigation - Restriction Enzyme Digestion** EDV0-KIT 225

In this experiment, the DNA from two suspects are each digested with two restriction enzymes in separate reactions and compared to crime scene samples after agarose gel electrophoresis. This flow chart outlines the procedure used for the restriction enzyme digestion of DNA obtained from Suspect 1. The DNA from Suspect 2 is digested in the same manner, using reaction tubes 3 and 4 (not shown).

**MODULE I OVERVIEW**

**Quick Reference:**

Dispensed Component	Tube Label
Crime scene DNA 1	CS 1
Crime scene DNA 2	CS 2
Suspect 1 DNA	DNA 1
Suspect 2 DNA	DNA 2
DNA Standard Marker	Markers
Enzyme Reaction Buffer	xxx buffer
Diluted Enzyme 1	Enzyme 1
Diluted Enzyme 2	Enzyme 2

Incubate at 45° C for 15 min. or at 37° C for 30-60 min.


Add 5 µl gel loading solution.

Samples are ready for electrophoresis.

To avoid cross-contamination, use a FRESH micropipet tip for each transfer of DNA and enzyme to the restriction enzyme reaction.

EDVOTEK® 1.800.EDVOTEK • Fax 202.370.1501 • info@edvotek.com • www.edvotek.com

Duplication of any part of this document is permitted for non-profit educational purposes only. Copyright © 2014 EDVOTEK, Inc. All rights reserved. 225 130913


# RESTRICTION ENZYMES

One of the most significant discoveries of molecular biology is a class of enzymes known as restriction endonucleases. These endonucleases (also known as restriction enzymes) are produced by many species of bacteria to protect themselves from invading viral DNA. Restriction enzymes act like molecular scissors, cutting double-stranded DNA at specific sequences. The utility of restriction enzymes has made molecular cloning, DNA mapping, sequencing and various genome-wide studies possible, launching the era of biotechnology.


Restriction Enzyme	Genus	Species	Strain	Recognition Site
<i>Ava I</i>	<i>Anabaena</i>	<i>variableis</i>	n/a	C <sup>*</sup> YGGG
<i>Bgl I</i>	<i>Bacillus</i>	<i>globigii</i>	n/a	GCCNMMN <sup>*</sup> NGGC
<i>EcoRI</i>	<i>Escherichia</i>	<i>coli</i>	RY 13	G <sup>*</sup> AATTC
<i>HaeIII</i>	<i>Haemophilus</i>	<i>aegyptius</i>	n/a	GG <sup>*</sup> CC
<i>HindIII</i>	<i>Haemophilus</i>	<i>influenzae</i>	R <sub>2</sub>	A <sup>*</sup> AGCTT
<i>Sac I</i>	<i>Streptomyces</i>	<i>achromogenes</i>	n/a	GAGCT <sup>*</sup> C

**Table 1:** Common Restriction Enzymes with Recognition Sites


Digestion by a restriction enzyme generates DNA fragments with one of two types of DNA ends—"sticky" or "blunt". To illustrate this, first consider the recognition site and cleavage pattern of *EcoRI*.


*EcoRI* cleaves between the G and neighboring A, as indicated by the arrows in the left side of the figure. It is important to note that the positions of the cleavage are staggered, so the resulting fragments project short overhangs of single-stranded DNA with complementary sequences. Such overhangs are referred to as "sticky" ends because the single-strands can interact with—or stick to—other overhangs with a complementary sequence (Figure 1). Digestion of the same piece of DNA using different enzymes can produce sticky ends of different lengths and strand orientation (5' vs. 3').

In contrast to *EcoRI*, *HaeIII* cuts both DNA strands at the same position, which generates fragments without an overhang. These so-called "blunt" ends can be joined with any other blunt end without regard for complementarity.

Some restriction enzymes, such as *AvaI*, recognize "degenerate" sites, which contain one or more variable positions.

Consequently, there are four possible sites that *AvaI* will recognize and cut: CCGGG, CCCGAG, CTCGGG and CTCGAG.

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